


DENTAL AMALGAM DISPOSAL

ONCE A WEEK: Inspect your amalgam separator and keep a log. Replace/Repair defects within 10 days of discovery.

DON'T USE BLEACH or chlorine containing cleaners to flush wastewater lines. Cleaners cannot have a pH lower than 6 or more than 8.



Keep amalgam out of sinks, toilets and trash.

Recycle chair-side traps, vacuum pump filters and amalgam separator waste.



Use precapsulated alloy not bulk.



Use ADA & ANSI compliant recycler.



City of Dallas



dallas water utilities
city of dallas

For more info, call 214-670-9725
DWUPretreatment@dallascityhall.com

© 2018 Goldstreet Design Agency, Inc. All Rights Reserved